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**\*\*English\*\***

**\*\*Warning:\*\*** Incorrect use of this device may cause burns, scars, pigmentation changes, or eye damage. Please read this user manual carefully before first use and keep it for future reference.

**\* \*\*Before Use:\*\*** Perform a skin test on a small area 48 hours before first use to check for adverse reactions. Do not use on tanned, sunburned, irritated, or skin with skin conditions (eczema, psoriasis, etc.). **\* \*\*During Use:\*\*** Wear the protective eyewear provided to protect your eyes. Do not use on tattoos, moles, warts, genital areas, or near the eyes. Avoid going over the same area multiple times in one session. **\* \*\*After Use:\*\*** Moisturize the treated skin with a soothing cream. Avoid sun exposure for at least 48 hours after treatment and apply sunscreen SPF 30 or higher. **\* \*\*General:\*\*** Keep this device out of reach of children. Do not use if pregnant or breastfeeding. If redness, irritation, or persistent pain occurs, discontinue use immediately and consult a doctor.

This document is automatically generated in response to European product safety regulations and in accordance with the GDPR.

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